


# Women's Health Intra-Uterine Contraceptive Devices

The Intra-Uterine contraceptive devices are ready for use in a sterile pack. They are formed of two flexible arms in polyethylene projecting from a vertical axis which is made radiopaque with barium sulphate.

## TT380<sup>®</sup> IUD

The TT380 IUD comes in two sizes, standard and short. Both devices are surrounded by a copper wire coil of 380mm<sup>2</sup> and fitted with a monofilament Nylon string at the base of the device.

Code: 7TT380  
Code: 7TT380Short



TT380

**NEW!**

## LOAD 375<sup>®</sup> IUD

The LOAD 375 IUD is standard size only. This device is surrounded by a copper wire coil of 375mm<sup>2</sup> and fitted with a monofilament Nylon string at the base of the device.

Code: 7LOAD375



LOAD 375

## Women's Health Products

### Bio-Cel Xcel Pregnancy Test

Rapid immunoassay for the qualitative detection of Human Chorionic Gonadotropin (hCG) in urine for the early detection of pregnancy.

- 2nd generation, one step pregnancy test
- Rapid: results in 3 minutes
- Easy to read colour
- Built-in test control
- Room temperature storage and long shelf life
- Self contained reagents - no mixing
- Single or batch testing - easy and rapid



### Emmett IUD Thread Retriever

Worldwide the Intra-Uterine Device (IUD) and the hormone-releasing Intra-Uterine System (IUS) are the most commonly used contraceptives.

At any one time 7-9% of IUD/IUS users report that the threads of the device can no longer be felt at the top of the vagina.

The Emmett IUD Thread Retriever is designed to retrieve the lost threads if they are within the endocervical canal or uterine cavity. The procedure is simple, painless, atraumatic – generally no analgesia nor local anaesthesia is required – and can be carried out in a Doctor's surgery or Clinic setting. This saves the patient time and inconvenience and may avoid the need for hysteroscopy.

Code: DUEMMETT, 5/PKT



DUEMMETT